

Dynamic Acquisition of ^{99m}Tc-Pyrophosphate SPECT/CT Images in Transthyretin Cardiac Amyloidosis: A Pilot Study

Journal of Nuclear Medicine 2025; 00:1–8 | DOI: 10.2967/jnumed.125.270405

Benjamin Auer, Alyssa De Moraes, Ardel J. Romero Pabon, Olivier F. Clerc, Sudhir Bhimaniya, Marie Foley Kijewski, Annu Kurian, Shilpa Vijayakumar, Sarah A.M. Cuddy, Marcelo F. Di Carli, and Sharmila Dorbala

Division of Nuclear Medicine and Molecular Imaging, Department of Radiology, Brigham and Women’s Hospital, Boston, Massachusetts

^{99m}Tc-pyrophosphate is commonly used to diagnose and stage Transthyretin Cardiac Amyloidosis (ATTR-CM). The study aimed to determine the timing of peak ^{99m}Tc-pyrophosphate uptake in the myocardium, blood pool, and bone, and to assess the feasibility of early imaging compared to standard late imaging in participants with and without ATTR-CM. Image analysis involved visual assessments and quantitative metrics, including SUVmean, SUVmax, and percentage injected dose per milliliter, alongside relative uptake ratios.






Figure 1. In the study design, the protocol takes advantage of VERITON-CT 4D (Dynamic SPECT/CT) quantitative capabilities as well as high sensitivity 3D imaging. Imaging includes a fast dynamic scan focused on the heart for the initial 5 minutes (for cardiac blood pool and perfusion), followed by a longer dynamic scan focused on the whole chest (for uptake kinetics and early washout kinetics). The 5-minute heart-focused dynamic scan achieved a temporal resolution of ~3 seconds. The subsequent 55-minute dynamic scan covered a larger FOV, resulting in a slightly coarser temporal resolution of ~5 seconds. Later phase quantitative SPECT/CT follows traditional imaging guidelines.



Blood pool, myocardium, background and bone regions of interest



4D and 3D SUV accuracy

Article Takeaways

- Suggests the possibility that early imaging may be possible when using quantitative SPECT/CT, improving the efficiency of ^{99m}Tc Pyrophosphate protocols for assessment of ATTR-CM
- Early quantitative ^{99m}Tc -PYP SPECT/CT imaging can reliably differentiate ATTR-CM from non-ATTR-CM cases, reducing the need for a prolonged 2-to-3-hour uptake period, streamlining clinical workflows, improving participant comfort, and supporting timely disease management.
- VERITON-CT provides accurate SUV in 4D and 3D SPECT/CT acquisitions. Quantitative dynamic SPECT/CT with VERITON-CT now enables kinetic evaluation of radiopharmaceutical accumulation in human clinical and research applications.




FIGURE 5. Comparison of absolute calibration-dependent and relative myocardial quantitative metrics between non-ATTR-CM and ATTR-CM groups for different periods. (A) Myocardial %ID/mL. (B) SUV_{mean} (left) and SUV_{max} (right) in $\text{LV}_{\text{myocardium}}$. (C) Semiquantitative SUV_{mean} ratios of myocardium to blood pool ($\text{LV}_{\text{myocardium}}/\text{LA}_{\text{top}}$, left) and myocardium to bone ($\text{LV}_{\text{myocardium}}/\text{LA}_{\text{top}}$, right). Data are presented as average \pm 95% CI, with individual values shown as dots. * $P \leq 0.05$. ** $P \leq 0.01$. *** $P \leq 0.001$. **** $P \leq 0.0001$. ns = not significant.



D-SPECT

CARDIO



VERITON-CT®



About Spectrum Dynamics Medical

Spectrum Dynamics Medical is spearheading the transformation of nuclear medicine from analog to digital technology, enabling clinicians to provide superior healthcare services with improved image quality and efficiency at lower doses.

D-SPECT®: the world's first digital cardiac SPECT system based on Broadview Technology: swiveling CZT digital detectors. It is the platform for advanced applications such as TruCorr Attenuation Correction, TruFlow 3D dynamic SPECT imaging for myocardial flow reserve (MFR) analysis and simultaneous dual isotope imaging (SDI).

VERITON-CT®: the world's first digital 360°CZT-based, ring-shaped gantry design SPECT/CT, with high-resolution 16 or 64 slice CT. Its 80cm wide bore accommodates patients of all sizes.

VERITON-CT continues to evolve with an energy range up to 400+ keV and multiple peak imaging, providing highest quality and quantitation-ready images.

Link to Article:

<https://jnm.snmjournals.org/content/early/2025/12/11/jnumed.125.270405>

Spectrum Dynamics Medical Information:

www.spectrum-dynamics.com

GlobalMarcom@spectrum-dynamics.com

Contact Us

Spectrum Dynamics Medical appreciates the opportunity to share our nuclear imaging solutions with you. To learn more about our products, services and how we can help you, please contact us.

Spectrum Dynamics
Medical HK Limited
Unit 1001, 10/F,
Mira Place Tower A
132 Nathan Road,
Tsim Sha Tsui
Kowloon, Hong Kong

Spectrum Dynamics
Medical SA
Rue de Lausanne, 31
1110 Morges Switzerland
Tel: +41 21 544 4710
Fax: +41 21 544 4711

Shanghai Guangmai
Medical Technology
Limited
708, Block A, No. 300,
Shuiyun Road, Pudong New
District, Shanghai, China.
Tel: 021-58281680

Spectrum Dynamics
Medical, Inc.
301 North Cattlemen
Road, Suite 109
Sarasota, FL 34232
Tel: +1 941 256 3541
Fax: +1 941 256 3832

Spectrum Dynamics
Medical Japan K.K.
Uchikanda TK building 6F
South, 1-5-13 Uchikanda
Chiyoda-ku, Tokyo Japan
Tel: +81 3 5843-9304
Fax: +81 3 5843-9305

